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
- DNA transfection
- siRNA transfection
- DNA/siRNA cotransfection
- CRISPR/Cas9
- Virus production

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
FACS analysis 24 h after GFP plasmid transfection according to the manufacturer's recommendation.

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
Phase contrast microscopy of HeLa cells 24 h after transfection according to the manufacturer's recommendation.

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
In vivo DNA & siRNA Delivery Reagent in vivo-jetPEI

IDEAL FOR SUBCUTANEOUS
TUMOR MODELS

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


Efficient and specific gene silencing following 1 nM INTERFERin-mediated siRNA transfection. Immunofluorescence microscopy 48 h after transfection

Silencing efficiency	
Adherent cell lines (1 nM siRNA)	
A549	> 90%
HeLa	
CaSki	
MCF7	
NIH-3T3	
RAW 264.7	
SiHa	60-70%
HepG2	
Primary cells (1 nM siRNA)	
Murine embryonic fibroblasts	> 90%
Primary human fibroblasts	
Primary human hepatocytes	
Suspension cell lines (5 nM siRNA)	
K562	> 80%
THP-1	

Selected successfully transfected cell lines and silencing efficiencies obtained with INTERFERin®.

Reagent	Reagent size	Buffer size	Polyplus Cat. No
in vitro			
jetPRIME Transfection reagent + buffer Broad range of adherent cell lines, DNA and/or siRNA transfection	0.1 ml	5 ml	114-01
	0.75 ml	60 ml	114-07
	1.5 ml	2 x 60 ml	114-15
	5 x 1.5 ml	10 x 60 ml	114-75
jetPRIME Transfection Reagent + buffer 5 x conc.	5 x 1.5 ml	120 ml	114-75C
INTERFERin Transfection Reagent siRNA day-to-day transfection	0.1 ml	-	409-01
	1 ml	-	409-10
	5 x 1 ml	-	409-50
FectoPRO™ DNA Transfection Reagent Bioproduction of Protein and Antibody	1 ml	1 ml	116-001
	10 ml	10 ml	116-010
	10 x 10 ml	10 x 10 ml	116-100
PEIpro™ DNA Transfection Reagent Protein & virus production	10 ml	-	115-010
	100 ml	-	115-100
	4 x 100 ml	-	115-400
jetPEI Transfection Reagent + NaCl solution DNA High Throughput Screening	0.1 ml	5 ml	101-01N
	1 ml	50 ml	101-10N
	4 x 1 ml	4 x 50 ml	101-40N
	10 ml	2 x 250 ml	101B-010N
jetPEI DNA Transfection Reagent DNA High Throughput Screening	1 ml	-	101-10
	4 x 1 ml	-	101-40
	10 ml	-	101B-010
jetPEI -Hepatocyte Transfection Reagent + NaCl solution DNA transfection into hepatocyte cells	0.1 ml	5 ml	102-01N
	0.5 ml	50 ml	102-05N
jetPEI -Macrophage Transfection Reagent + NaCl solution DNA transfection into macrophage cells	0.1 ml	5 ml	103-01N
	0.5 ml	50 ml	103-05N
jetPEI -HUVEC Transfection Reagent + NaCl solution DNA transfection into HUVEC	0.1 ml	5 ml	108-01N
	0.5 ml	50 ml	108-05N
PULSin™ Protein, Antibody and Peptide Delivery Reagent Hepes solution + 200 µL positive control included	0.1 ml	20 ml	501-01
	0.4 ml	20 ml	501-04
INTERFERin -HTS Transfection Reagent siRNA High Throughput Screening	0.2 ml	-	410-002
	1.5 ml	-	410-015
in vivo			
in vivo-jetPEI Delivery Reagent + glucose sol. DNA and siRNA in vivo delivery	0.1 ml	10 ml	201-10G
	1.5 ml	2 x 10 ml	201-50G
in vivo-jetPEI -Gal Delivery Reagent + glucose sol., Galactose-conjugated	0.1 ml	10 ml	202-10G
in vivo-jetPEI -Man Delivery Reagent + glucose sol., Mannose-conjugated	0.1 ml	10 ml	203-10G
jetSI™ 10 mM Delivery Reagent, siRNA delivery into the brain	0.5 ml	-	403-05

Polyplus-transfection started in 2001 to market transfection products developed in the lab of Dr. Jean-Paul Behr, who is considered one of the preminent transfection researchers in the world. Seminal work was published in 1989 on the use of polycationic lipids for the use of transfection⁽¹⁾ that led to several still popular transfection reagents. The next advancement in this field is associated with Polyethylenimine (PEI), a cationic polymer. PEI was discovered as an improved transfection reagent in 1995⁽²⁾ and its associated «proton sponge» transfection paradigm is at the root of many modern transfection reagents.

Polyplus-transfection continues its strong tradition of scientific excellence and innovation offering the best reagents for delivery of nucleic acids in the market today.

(1) J.-P. Behr, B. Demeneix, J.-P. Loeffler, J. Perez-Mutul, «Efficient gene transfer into mammalian primary endocrine cells with lipopolyamine coated DNA», Proc. Natl. Acad. Sci. USA, 86, 6982, 1989.
(2) O. Boussif, P. Lezoualc'h, M.A. Zanta, M.D. Mergny, D. Scherman, B. Demeneix, J.P. Behr, «A versatile vector for gene and oligonucleotide transfer into cells in culture and in vivo: polyethylenimine.» Proc. Natl. Acad. Sci. USA, 92, 7297-7301, 1995.

Prices and availability are subject to change. Copyright© 2015 Polyplus-transfection SA. All Rights Reserved. Each product is sold with a limited warranty which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. Polyplus-transfection SA products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Polyplus-transfection SA. INTERFERin, jetPEI, jetPRIME, jetSI, PEIpro, Polyplus-transfection and PULSin are trademarks or registered trademarks of Polyplus-transfection SA. Opti-MEM is a registered trademark of Life Technologies.