

# *A little bit of everything* for cell research




*Just capture images, make movies*

*as well as check the expression level*


## Real time cell growth curve

HeLa cells growth were observed for 40 hours with 10 minutes intervals and analyzed monolayer confluence. For apoptosis assays, experimental group (channel 2) was treated with Sataurosporine.




## Cell migration (wound healing) assay

U-2 OS stable cells were incubated for 100 hours after scratch. JuLI™ FL calculated the confluence with growth of surface unfarmed automatically.



## Fluorescence expression level checking

JuLI™ FL counting software provides actual cell counts, viability and automatically calculates cell concentration. It also can identify and measure the levels of fluorescence (GFP or RFP expression).



# Meet JuLI™ FL to meet your needs for live cell imaging

The JuLI™ FL, the fluorescence cell history recorder, with the states-of-the-art optics enables users to perform the image capturing and the time-lapse movie making. It also can support cell counting and measuring the level of fluorescence expression.

## Automated cell confluence detection

Quantified cell confluence results with low variation



## 10.1" color LCD touch screen

User friendly interface




## Dual system (\*Optional)

Compare control and experimental samples using dual system, concurrently




## Live cell image capture

U-2 OS cells were transfected using the Neon™ transfection system (from Life Technologies) and 0.25 µg of the EGFP-N1 plasmid & RFP plasmid. 48 hours after transfection, the images are captured.



## Time-lapse image capture and movie making

Cell-growth images were captured for 63 hours with 20 minutes intervals in U-2 OS GFP & RFP stable cell line.






# It works inside incubator!

## Applications

- Cell growth monitoring
- Cell migration assay
- Cell confluence detection
- Cell viability & counting
- Cell culture quality control
- Proliferation assay
- Cell-based assay optimization



Station unit



Scope unit

## Specifications

Magnification	4 X and digital zoom ( ~ 450 X)
Image resolution	1280 x 960 pixels (1.3 M)
Exported formats	JPEG, TIFF, BMP, PNG (image), AVI (movie), CSV (raw data)
Display	10.1" LCD touch screen
Light source	GFP channel (JULI-FLG04): Blue LED RFP channel (JULI-FLR04): Green LED
Optical filter	GFP channel (JULI-FLG04): Excitation 466/40 nm, Emission 525/50 nm RFP channel (JULI-FLR04): Excitation 525/50 nm, Emission 580LP
Dimension & Weight	Scope unit : 300 x 190 x 188 mm, 4.5 kg Station unit : 282 x 285 x 160 mm, 3.2 kg
Storage	320 GB Hard drive, 4 GB USB drive
Power	AC100 - 240Va.c., 1.9A, 100w, 50/60 Hz
Operating Environment	5 - 40 °C, 20 - 95 %

## Ordering information

	Cat. No.	Description	Remark
Device	JULI-FLG04	JuLI FL (Single set)	For GFP channel
	JULI-FLGD04	JuLI FL (Dual set)	
	JULI-FGSC	Scope unit of JuLI FL	
	JULI-FLR04	JuLI FL (Single set)	For RFP channel
	JULI-FLRD04	JuLI FL (Dual set)	
	JULI-FRSC	Scope unit of JuLI FL	
Accessories	JULI-BRCM	Counting starter kit	Counting module set, Cell counting slide
	JULI-BRTB	XY stage	
Disposable	JULI-BRS50	Cell counting slide	50 slides (100 counts)
	JULI-BRS1000	Cell counting slide	1,000 slides (2,000 counts)
	JULI-BRS5000	Cell counting slide	5,000 slides (10,000 counts)

NanoEnTek

NESCT-JUF-001E(V.1.0)

Demo request

demo@nanoentek.com

website

www.nanoentek.com

e-mail

sales@nanoentek.com

### NanoEnTek Inc.

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul,  
08389, Korea  
Tel : +82-2-6220-7940  
Fax : +82-2-6220-7721

### NanoEnTek USA Inc.

5627 Stoneridge Drive Suite 304, Pleasanton,  
CA 94588, USA  
Tel : +1-925-225-0108, +1-888-988-0108(Toll free)  
Fax : +1-925-225-0109